

# **SYSTEMATIC ANALYSIS OF INORGANIC SALT MIXTURE**

**II B.Sc. CHEMISTRY PRACTICALS**

**PG DEPARTMENT OF CHEMISTRY  
SWAMY DAYANANDA COLLEGE OF ARTS &  
SCIENCE, MANJAKKUDI-612610, Tiruvarur-Dt**

By Dr.G.Marimuthu

## SYSTEMATIC INORGANIC QUALITATIVE ANALYSIS

The given unknown inorganic mixture contains two inorganic salts. Each salt contains anion and a cation. One of the two anions, is an interfering anion and the other one is a non-interfering anion.

The list of non-interfering anion is given below

1. Carbonate anion      2. Sulphide anion      3. Sulphate anion      4. Nitrate anion      5. Chloride anion      6. Bromide anion

The list of interfering anions, is given below

1. Fluoride anion      2. Borate anion      3. oxalate anion      4. Phosphate anion

The group wise list of cations, is given below :

Group I	Group II	Group III	Group IV	Group V	Group VI
1. Lead	1. Copper 2. Bismuth 3. Cadmium	1. Iron 2. Aluminium	1. Cobalt 2. Nickel 3. Manganese 4. Zinc	1. Barium 2. Strontium 3. Calcium	1. Ammonium 2. Magnesium

Analysis of anions must be done first followed by the analysis of cations

It is necessary to do all the available tests for the anions

It is essential to record all the alterations as and when they are made

Negative reactions should not be ignored and should be recorded

### PRELIMINARY TESTS

#### 1. Visual Examination of the Sample

NO	EXPERIMENT	OBSERVATION	INFERENCE
I)	COLOUR The colour of the given mixture is examined	a) Black b) Blue c) Brown d) Green e) Pink f) No characteristic colour	May be Pb, Cu, Fe or Ni salt May be Cu or Fe salt May be Pb, Fe or Bi salt May be Fe, Cu or Ni salt May be Co salt May not be Pb, Bi, Cu, Fe, C, Ni etc.

I) i)	<p><b>ACTION OF HEAT</b></p> <p>A small amount of the given mixutre, is taken in a dry test tube and heated strongly</p>	<p>a) Becomes black</p> <p>b) Yellow when hot and white when cold</p> <p>c) Colourless and odourless gas, turning lime water milky  <math>\text{CO}_3^{2-} \rightarrow \text{O}_2^{2-} + \text{CO}_2</math></p> <p><math>\text{CO}_2 + \text{Ca(OH)}_2 \rightarrow \text{H}_2\text{O} + \text{CaCO}_3</math> (milky)</p> <p>d) Colourless gas with the smell of burnt sulphur, turning potassium dichromate paper, green  <math>2\text{S}^{2-} + 3\text{O}_2 \rightarrow 2\text{O}^{2-} + 2\text{SO}_2</math> (smell of burnt S]  <math>3\text{SO}_2 + \text{Cr}_2\text{O}_7^{2-} + 8\text{H}^+ \rightarrow 3\text{SO}_3 + 4\text{H}_2\text{O} + 2\text{Cr}^{3+}</math> (green)</p> <p>e) Colourless gas with the smell of burnt sulphur, turning potassium dichromate paper, green  <math>\text{SO}_3^{2-} \rightarrow \text{O}^{2-} + \text{SO}_2</math> (smell of burnt sulphur),  <math>3\text{SO}_2 + \text{Cr}_2\text{O}_7^{2-} + 8\text{H}^+ \rightarrow 3\text{SO}_3 + 4\text{H}_2\text{O} + 2\text{Cr}^{3+}</math> (green)</p> <p>f) Colorless gas with the smell of ammonia  <math>\text{NH}_4^+ \rightarrow \text{H}^+ + \text{NH}_3</math> (smell of ammonia)</p> <p>g) Colourless gas, burning with a blue flame  <math>(\text{COO})_2^{2-} \rightarrow \text{O}^{2-} + \text{CO} + \text{CO}_2</math> (blue flame)</p> <p>h) Reddish-brown gas, turning ferrous sulphate paper, brown  <math>4\text{NO}_3^- \rightarrow 2\text{O}_2 + \text{O}_2 + 4\text{NO}_2</math> (reddish-brown)</p> <p>(I) No characteristic reaction</p>	<p>May be copper salt  May be zinc salt  May be carbonate or oxalate anion</p> <p>May be sulphide anion</p> <p>May be sulphate anion</p> <p>May be ammonium salt</p> <p>May be oxalate anion</p> <p>My be nitrate anion</p> <p>May not be carbonate, sulphide, sulphite, oxalate nitrate etc.</p>
----------	--	--	---

ii)	<p><b>FLAME COLOURATION</b></p> <p>a) a small amount of the given mixture is mixed with a drop of conc. HCl to make a paste. This paste is introduced at the base of the outer oxidizing region of the non-luminous flame of the bunsen burner</p>	<p>a) Brick – red flame  b) Crimson-red flame  c) Green flame  d) No characteristic colouration</p>	<p>May be calcium salt  May be strontium salt  May be Ba, Cu or borate salt  Absence of Ba, Ca and Sr salt</p>
iii)	<p>A small amount of the given mixture is mixed with an equal amount of calcium fluoride salt and a drop of Conc H<sub>2</sub>SO<sub>4</sub> to make a paste. This paste is introduced at the base of the outer oxidizing region of the non-luminous flame of the bunsen burner.</p>	<p>Green-edged flame appears  <math>\text{CaF}_2 + \text{H}_2\text{SO}_4 \rightarrow \text{CaSO}_4 + 2 \text{HF}</math>  <math>\text{B}_2\text{O}_3 + 3 \text{HF} \rightarrow 3 \text{H}_2\text{O} + 2\text{BF}_3</math>  (green flame)</p>	<p>Presence of borate anion</p>
iv)	<p>A small amount of the given mixture is mixed with an equal amount of borax salt and a drop of conc. H<sub>2</sub>SO<sub>4</sub> to make a paste. This paste is introduced at the base of the outer oxidizing region of the non-luminous flame of the bunsen burner.</p>	<p>Green edged flame appears</p> <p>No green-edged flame appears</p>	<p>Presence of fluoride anion</p> <p>Absence of fluoride anion</p>
v)	<p><b>Ethyl borate test</b></p> <p>A small amount of the given mixture is taken in a test tube and heated with a few drops of ethanol and conc H<sub>2</sub>SO<sub>4</sub>. the gas evolved, is introduced at the base of the outer oxidizing region of the non-luminous flame of the bunsen burner.</p>	<p>Green edged flame appears  <math>\text{BO}_3^{3-} + 3\text{H}^+ \rightarrow \text{H}_3\text{BO}_3</math>  <math>\text{H}_3\text{BO}_3 + 3 \text{EtOH} \rightarrow (\text{Et})_3\text{BO}_3 + 3 \text{H}_2\text{O}</math></p> <p>No green-edged flame</p>	<p>Presence of borate anion</p> <p>Absence of borate anion</p>
2.	<p><b>PRELIMINARY EXAMINATION FOR ANIONS</b></p> <p><b>ACTION OF DILUTE HCl</b></p> <p>i) A small amount of the given mixture is taken in a test tube and treated / heated with a few drops of dilute HCl</p>	<p>(a) Brisk effervescence of a colourless, odourless gas, turning lime water, milky.  <math>\text{CO}_3^{2-} + 2\text{H}^+ \rightarrow \text{H}_2\text{O} + \text{CO}_2</math>  (Colourless gas)  <math>\text{CO}_2 + \text{Ca}(\text{OH})_2 \rightarrow \text{H}_2\text{O} + \text{CaCO}_3</math>  (milky)</p> <p>b) Colourless gas with rotten-egg smell, turning lead acetate paper,</p>	<p>May be carbonate anion</p> <p>May be sulphide anion</p>

		<p>black and alkaline nitroprusside paper red.</p> $\text{S}^{2-} + 2\text{H}^+ \rightarrow \text{H}_2\text{S}$ $\text{H}_2\text{S} + \text{Pb}(\text{OAc})_2 \rightarrow \text{HOAc} + \text{Pb S}$ <p>(black)</p>	
		<p>c) Colourless gas with the smell of burnt sulphur, turning acidified potassium dichromate paper, green</p> $\text{SO}_3^{2-} + 2\text{H}^+ \rightarrow \text{H}_2\text{O} + \text{SO}_2 \text{ (smell of burnt S)}$ $3\text{SO}_2 + \text{CrO}_7^{2-} + 8\text{H}^+ \rightarrow 3\text{SO}_3 + 4\text{H}_2\text{O} + 2\text{Cr}_3^+ \text{ (green)}$ <p>d) Reddish brown gas evolves and the solution becomes blue in colour</p> $\text{NO}_2^- + \text{H}^+ \rightarrow \text{HNO}_2$ $3\text{HNO}_2 \rightarrow \text{H}_2\text{O} + \text{HNO}_3 + 2\text{NO}$ $2\text{NO} + \text{O}_2 \rightarrow 2\text{NO}_2 \text{ (reddish-brown gas)}$ $\text{NO}_3^- + \text{H}^+ \rightarrow \text{HNO}_3$ $4\text{HNO}_3 \rightarrow \text{H}_2\text{O} + \text{O}_2 + 4\text{NO}_2 \text{ (reddish-brown gas)}$ <p>d) No characteristic reaction</p>	<p>May not be carbonate, sulphide and sulphite anion</p>
ii)	ACTION OF CON. HCl	<p>a) Brisk effervescence of a colourless, odourless gas, turning lime water, milky.</p> $\text{CO}_3^{2-} + 2\text{H}^+ \rightarrow \text{H}_2\text{O} + \text{CO}_2 \text{ (Colourless gas)}$ $\text{CO}_2 + \text{Ca}(\text{OH})_2 \rightarrow \text{H}_2\text{O} + \text{Ca CO}_3 \text{ (milky)}$ <p>b) Colourless gas with rotten-egg smell, turning lead acetate paper, black and alkaline nitroprusside paper red.</p> $\text{S}^{2-} + 2\text{H}^+ \rightarrow \text{H}_2\text{S}$ $\text{H}_2\text{S} + \text{Pb}(\text{OAc})_2 \rightarrow \text{HOAc} + \text{Pb S}$ <p>(black)</p>	<p>May be carbonate anion</p> <p>May be sulphide anion</p>

		<p>c) Colourless gas with the smell of burnt sulphur, turning acidified potassium dichromate paper, green</p> $\text{SO}_3^{2-} + 2\text{H}^+ \rightarrow \text{H}_2\text{O} + \text{SO}_2 \text{ (smell of burnt S)}$ $3\text{SO}_2 + \text{CrO}_7^{2-} + 8\text{H}^+ \rightarrow 3\text{SO}_3 + 4\text{H}_2\text{O} + 2\text{Cr}_3^+ \text{ (green)}$ <p>d) No characteristics reaction</p>	May not be carbonate, sulphide and sulphite anion
iii)	<p><b>ACTION OF DILUTE H<sub>2</sub>SO<sub>4</sub> AND KMnO<sub>2</sub></b> A small amount of the given mixture is taken in a test tube and treated with a few drops of dilute H<sub>2</sub>SO<sub>4</sub> and KMnO<sub>4</sub> solution.</p>	<p>(a) Pink colour disappears immediately in cold condition</p> $2\text{MnO}_4^- + 5\text{NO}_2^- + 8\text{H}^+ \rightarrow 5\text{NO}_3^- + 3\text{H}_2\text{O} + 2\text{Mn}^{2+}$ <p>(b) Pink colour disappears only slowly.</p> <p>(c) Pink colour disappears only on heating</p> $2\text{MnO}_4^- + 16\text{H}^+ + 5(\text{COO})_2^{2-} \rightarrow 10\text{CO}_2 + 8\text{H}_2\text{O} + 2\text{Mn}^{2+}$ <p>d) No characteristics reaction</p>	May not be carbonate, sulphide and sulphite anion
iv)	<p><b>ACTION OF DILUTE H<sub>2</sub>SO<sub>4</sub> AND MnO<sub>4</sub></b> A small amount of the given mixture is mixed with an equal amount of MnO<sub>2</sub> powder and heated with dilute H<sub>2</sub>SO<sub>4</sub> in a test tube.</p>	<p>(a) Brisk effervescence of a colourless, odourless gas, turning lime water, milky</p> $\text{CO}_3^{2-} + 2\text{H}^+ \rightarrow \text{H}_2\text{O} + \text{CO}_2 \text{ (colourless gas)}$ $\text{CO}_2 + \text{Ca(OH)}_2 \rightarrow \text{H}_2\text{O} + \text{CaCO}_3 \text{ (milky)}$ <p>b) Brisk effervescence of a colourless, odourless gas, turning lime water, milky</p> $\text{MnO}_2 + 4\text{H}^+ + (\text{COO})_2^{2-} \rightarrow \text{Mn}^{2+} + 2\text{H}_2\text{O} + 2\text{CO}_2$ $\text{CO}_2 + \text{Ca(OH)}_2 \rightarrow \text{H}_2\text{O} + \text{CaCO}_3 \text{ (milky)}$ <p>(c) No characteristic reaction</p>	<p>Presence of carbonate anion</p> <p>May be oxalate anion (colourless gas)</p> <p>May not be carbonate or oxalate anion</p>
V	<p><b>ACTION OF CONC H<sub>2</sub>SO<sub>4</sub> AND MnO<sub>2</sub></b> A small amount of the given mixture is mixed with an equal amount of MnO<sub>2</sub> powder and heated with conc. H<sub>2</sub>SO<sub>4</sub> in a test tube</p>	<p>(a) Greenish-yellow gas with pungent smell, turning starch potassium iodide paper, blue</p>	<p>May be chloride anion</p> <p>May be bromide anion</p>





	<p>A small amount of the given mixture is dissolved in dilute H<sub>2</sub>SO<sub>4</sub> in a test tube and a few drops of dilute AgNO<sub>3</sub> solution are added.</p>	<p>a) A white precipitate, soluble in dilute HNO<sub>3</sub> forms a black precipitate on boiling  <math>\text{SO}_3^{2-} + 2 \text{Ag}^+ \rightarrow \text{Ag}_2 \text{SO}_3</math> (white precipitate)  2. <math>\text{Ag}_2 \text{SO}_3 \rightarrow \text{O}_2 + 2 \text{SO}_2 + 4 \text{Ag}</math> (black precipitate)  b) A black precipitate, insoluble in cold dilute HNO<sub>3</sub> but dissolves on heating.  <math>\text{S}^{2-} + \text{AgNO}_3 \rightarrow \text{NO}_3^- + \text{Ag}_2\text{S}</math> (black precipitate)  (c) No characteristic reaction</p>	<p>May be sulphite anion</p> <p>May be sulphide anion</p> <p>May not be sulphide or sulphite anion</p>
ix)	<p><b>ACTION OF BaCl<sub>2</sub> SOLUTION</b></p> <p>A small amount of the given mixture is dissolved in dilute HCl in a test tube and a few drops of dilute BaCl<sub>2</sub> solution are added</p>	<p>a) A white precipitate soluble in dilute HCl  <math>\text{SO}_3^{2-} + \text{BaCl}_2 \rightarrow 2 \text{Cl}^- + \text{BaSO}_3</math> (white precipitate)  b) A white precipitate, insoluble in hot dilute HCl and dilute HNO<sub>3</sub>  <math>\text{SO}_4^{2-} + \text{BaCl}_2 \rightarrow 2 \text{Cl}^- + \text{BaSO}_4</math> (white precipitate)  (c) No characteristic reaction</p>	<p>May be sulphite anion</p> <p>May be sulphate anion</p> <p>May not be sulphite or sulphate anion</p>
x)	<p><b>CHROMYL CHLORIDE TEST</b></p> <p>A small amount of the given mixture is taken in a test tube and heated strongly with a little of Potassium dichromate and conc. H<sub>2</sub>SO<sub>4</sub></p>	<p>a) Orange-red vapours. condensing to a red liquid forms yellow precipitate with Pb (OAc)<sub>2</sub></p> $\text{Cl}^- + \text{H}^+ \rightarrow \text{HCl}$ $\text{Cr}_2\text{O}_7^{2-} + 2\text{H}^+ \rightarrow \text{H}_2\text{O} + 2 \text{CrO}_3$ $\text{CrO}_3 + 2\text{HCl} \rightarrow \text{H}_2\text{O} + \text{CrO}_2\text{Cl}_2$ $\text{CrO}_2\text{Cl}_2 + 2\text{H}_2\text{O} \rightarrow 2 \text{HCl} + \text{H}_2\text{CrO}_4$ $\text{H}_2\text{CrO}_4 + \text{Pb} (\text{OAc})_2 \rightarrow 2 \text{HOAc} + \text{Pb CrO}_4 \text{ (Yellow precipitate)}$	<p>Presence of chloride anion</p> <p>Absence of chloride anion</p>

		(b) No characteristic reaction	
xi)	<b>ACTION OF CONC HNO<sub>3</sub></b> A small amount of the given mixture is taken in a test tube and warmed with a few drops of conc. HNO <sub>3</sub> till no more brown fumes evolve and then 2 cc of ammonium molybdate reagent is added to it.	a) Canary yellow precipitate appears in cold on shaking or on slight heating on vigorous boiling b) No characteristic reaction	Presence of phosphate anion  Absence of phosphate or arsenate anions
xii)	<b>ACTION OF NEUTRAL FeCl<sub>3</sub></b> A small amount of the given mixture is taken in a test tube and warmed with a few drops of dilute HCl and a few drops of Fe Cl <sub>3</sub> are added	a) A pale – yellow precipitate appears $\text{PO}_4^{3-} + \text{FeCl}_3 \rightarrow 3 \text{Cl}^- + \text{FePO}$ (pale-yellow precipitate) b) No characteristic reaction	Presence of phosphate anion  Absence of phosphate anion
<b>3. EXAMINATION OF ANIONS IN SOLUTION PREPARATION OF SODIUM CARBONATE EXTRACT</b>			
	A small amount of the given mixture and sodium carbonate salt are taken in the ratio of 1:3 in a clean conical flask and a test tube full of distilled water is added to it. Then the content is boiled for 5 minutes and filtered in a test tube to get the extract. The extract is divided into several portions and the following tests are carried out.		
(I)	<b>PRECIPITATION OF SILVER SALT</b> A portion of the extract is neutralized with dilute HNO <sub>3</sub> added dropwise with shaking till the effervescence ceases and centrifuged. To the centrifugate, a few drops of silver nitrate reagent are added and centrifuged.	a) A curdy-white precipitate appears, dissolves readily in aqueous ammonia and reprecipitates with dilute HNO <sub>3</sub> $\text{Cl}^- + \text{AgNO}_3 \rightarrow \text{NO}_3^- + \text{AgCl}$ (curdy-white precipitate) b) A pale-yellow precipitate appears, dissolves readily in aqueous ammonia and reprecipitates with dilute HNO <sub>3</sub> $\text{Br}^- + \text{AgNO}_3 \rightarrow \text{NO}_3^- + \text{Ag Br}$ (Pale-yellow precipitate)	May be chloride anion  Presence of Bromide Anion
	To the centrifugate, aqueous Ammonia is added gradually, along the sides of the test tube without shaking.	a) Yellow ring appears at the junction of solutions, soluble in NH <sub>4</sub> OH and in dilute HNO <sub>3</sub> $\text{HNO}_3 \rightarrow \text{PO}_4^{3-} + 3 \text{AgNO}_3 \quad 3 \text{NO}_3^- + \text{Ag}_3 \text{PO}_4$ (Yellow ring) b) No characteristic reaction	Presence of phosphate  Absence of chloride, bromide, and phosphate anions

ii)	<p><b>PRECIPITATION OF BARIUM SALTS</b></p> <p>A portion of the extract is neutralized with dilute HCl, added dropwise with shaking till the effervescence ceases and centrifuged. To the centrifugate, a few drops of Barium chloride are added and centrifuged.</p>	<p>a) White precipitate appears, insoluble in hot dilute HCl and dilute HNO<sub>3</sub>  <math>\text{SO}_4^{2-} + \text{BaCl}_2 \rightarrow 2 \text{Cl}^- + \text{BaSO}_4</math> (white precipitate)</p> <p>b) No white precipitate appears</p>	<p>Presence sulphate anion</p> <p>Absence of sulphate Anion</p>
	<p>To the centrifugate, aqueous ammonia is added and centrifuged again</p>	<p>a) White precipitate appears and dissolves in dilute acetic acid</p> <p>b) While precipitate appears and found insoluble in dilute acetic acid</p> <p>c) No characteristic reaction</p>	<p>Presence of phosphate or borate anion</p> <p>Presence of oxalate or fluoride anion</p> <p>Absence of phosphate, borate, oxalate and fluoride anion</p>
iii)	<p><b>PRECIPITATION OF LEAD SALT</b></p> <p>A portion of the extract is neutralized with dilute acetic acid, added dropwise with shaking till the effervescence ceases and centrifuged. To the centrifugate, a few dorps of lead and acetate solution are added and centrifuged.</p>	<p>a) A white precipitate appears, dissolves readily in ammonium acetate anion  <math>\text{SO}_3^{2-} + \text{Pb}(\text{OAc})_2 \rightarrow 2 \text{OAc}^- + \text{PbSO}_4</math> (white precipitaate)</p> <p>b) No characteristic reaction</p>	<p>Presence of sulphate anion</p> <p>Absence of sulphate anion</p>
iv)	<p><b>PRECIPTATION OF CALCIUM SALTS</b></p> <p>A portion of the extract is neutralized with dilute acetic acid added drop wise with shaking till the effervescence ceases and centrifuged. To the centrifuge a few drops of calcium chloride solution are added and centrifuged.</p> <p>To the residue 2 drops of dilute H<sub>2</sub>SO<sub>4</sub> and a dorp of very dilute kmNO<sub>4</sub> are added.</p>	<p>a) White precipitate, does not decolourise pink colour of KMnO<sub>4</sub> solution, soluble in HOAc  <math>\text{F}^- + \text{CaCl}_2 \rightarrow 2 \text{Cl}^- + \text{CaF}_2</math> (white precipitate)</p> <p>b) White precipitate, decolourises pink colour of KMnO<sub>4</sub> solution  <math>(\text{COO})_2^{2-} + \text{CaCl}_2 \rightarrow 2 \text{Cl}^- + \text{Ca}(\text{COO})_2</math> (white precipitate)</p> <p>c) No characteristic reaction</p>	<p>Presence of fluoride anion</p> <p>Presence of oxalate anion</p> <p>Absence of fluoride and oxalate anions</p>
v)	<p><b>BROWN RING TEST</b></p>		<p>Presence of nitrate anion</p>

	<p>A portion of the extract neutralized with dilute H<sub>2</sub>SO<sub>4</sub> added dropwise with shaking till the effervescence ceases and a freshly prepared solution of ferrous sulphate is added without shaking. Then a drop of conc H<sub>2</sub>SO<sub>4</sub> is added along the side of the test tube</p>	<p>a) A brown ring appears at the junction of the 2 layers  <math>\text{NO}_3^- + \text{H}^+ \rightarrow \text{HNO}_3</math>  <math>\text{Fe}^{2+} + \text{NO}_3^- + 4\text{H}^+ \rightarrow \text{Fe}^{3+} + 2\text{H}_2\text{O} + \text{NO}</math>  <math>\text{FeSO}_4 + \text{NO} + 5\text{H}_2\text{O} \rightarrow [\text{Fe}(\text{H}_2\text{O})_3\text{NO}]\text{SO}_4</math> (brown ring)          b) No brown ring appears</p>	<p>Absence of nitrate anion</p>
<p><b>SYSTEMATIC ANALYSIS FOR CATIONS:</b>  <b>IN SOLUTION ELIMINATION OF INTERFERING ANIONS</b></p> <p>Interfering anions are anions which cause precipitation of cations outside their normal groups. Salts of such anions are soluble in mineral acids, but are insoluble in neutral or alkaline solution. Since dilute HCl is the group reagent evolved for the analysis of groups I and II, interfering anions will not affect the analysis of group I and II cations.</p> <p>Since aqueous ammonia (base) is the group reagent involved for the analysis of group III, interfering anions will cause the precipitation of not only group III cations, but also the cations of group IV, V and VI.</p> <p>So, it is essential to eliminate the interfering anions, before entering into the systematic analysis of cations</p> <p><b>ELIMINATION OF INTERFERING OXALATE ANION</b></p> <p>A little of the given unknown mixture is roasted for 15 minutes in a china dish. The residue is used for group separation analysis.</p> <p><b>ELIMINATION OF INTERFERING FLUORIDE OR BORATE ANION</b></p> <p>A little of the given unknown mixture is heated in a china dish with 1 cc of conc. HCl to dryness. The process is repeated 4 to 5 times. The final residue is utilized for group separation analysis.</p> $\text{F}^- + \text{H}^+ \rightarrow \text{HF}$ $\text{BO}_3^- + 3\text{H}^+ \rightarrow \text{HBO}_3$ <p><b>ELIMINATION OF INTERFERING PHOSPHATE ANION</b></p> <p>Elimination of phosphate anion is to be done at the beginning of group III analysis. The residue for group III contains actually phosphates of non-group III cations as well. The centrifugate from group II is boiled well, to expel all hydrogen sulphide and a drop of conc. HNO<sub>3</sub> acid is added and boiled gently. Then a few drops of ammonium chloride solution and 1 ml of zirconyl chloride reagent are added and digested for 2 minutes it is centrifuged and the centrifugate is tested for the presence of phosphate anion by adding zirconyl chloride reagent again. The same process is repeated till all phosphate has been precipitated and this residue of zirconium phosphate is rejected. The centrifugate is used for group separation analysis from group III.</p>			

A small amount of the given unknown mixture is made to dissolve in a suitable solvent like dilute or conc. HCl or H <sub>2</sub> SO <sub>4</sub> or HNO <sub>3</sub> . to this solution, group I reagent, namely dilute HCl is added to slight excess, for complete precipitation and centrifuged. $Pb^{2+} + 2 HCl \rightarrow 2 H^+ + PbCl_2$ group I precipitate)			
RESIDUE Presence of group I cation PbCl <sub>2</sub>	CENTRIFUGATE for group II analysis To the centrifugate, group II reagent namely hydrogen sulphide gas passed to saturation and centrifuged. $Cu^{2+} + H_2S \rightarrow 2 H^+ + Cu S$ (black precipitate for group II) $Bi^{3+} + 3 H_2S \rightarrow 6H^+ + Bi_2 S_3$ (precipitate for group II) $Cd^{2+} + H_2S \rightarrow 2H^+ + Cd S$ (precipitate for group II)		
NO RESIDUE Absence of Group I cation	RESIDUE Presence of group II cations Cu S or Bi <sub>2</sub> S <sub>3</sub> or Cd S NO RESIDUE Absence of Group II cations	CENTRIFUGATE (for group III analysis) The centrifugate is boiled to expel all excess of hydrogen sulphide gas and a drop of HNO <sub>3</sub> is added and boiled for sometime. Then group III reagents namely ammonium chloride solution and aqueous ammonia (slight excess) are added and centrifuged. $Mn^{2+} + 2 NH_4 OH \rightarrow 2 NH_4^+ + Mn (OH)_2$ $Co^{2+} + H_2S \rightarrow 2 H^+ + CoS$ $Zn^{2+} + H_2S \rightarrow 2 H^+ + Zn S$ (precipitate for group IV) $Ni^{2+} + H_2S \rightarrow 2 H^+ + Ni S$	
		RESIDUE Presence of group IV cations Mn S, CO S, Zn S and Ni S NO RESIDUE Absence of group IV cations	CENTRIFUGATE (for group V analysis) The centrifugate is boiled to expel all excess of hydrogen sulphide gas and a few drops of dilute HNO <sub>3</sub> are added and again boiled to concentrate the solution. Then group v reagent namely ammonium carbonate solution is added for complete precipitation and centrifugate. $Ca^{2+} + (NH_4)_2CO_3 \rightarrow 2 NH_4^+ + CaCO_3$ $Ba^{2+} + (NH_4)_2CO_3 \rightarrow 2 NH_4^+ + Ba CO_3$ (group IV precipitate) $Sr^{2+} + (NH_4)_2CO_3 \rightarrow 2 NH_4^+ + SrCO_3$
			RESIDUE Presence of group V cations Ca CO <sub>3</sub> , Ba CO <sub>3</sub> and Sr CO <sub>3</sub> NO RESIDUE
			CENTRIFUGATE (for group VI analysis) The centrifugate is used for the analysis of group VI cations

			Absence of group V cation	
--	--	--	------------------------------	--

### ANALYSIS OF GROUP I CATIONS

The residue for group I is washed with water and boiled for sometime and centrifuged	
<b>NO RESIDUE</b> Absence of Mercury and silver cations	<b>CENTRIFUGATE :</b> 1. To one portion of the centrifugate, potassium chromate solutions is added <div style="display: flex; justify-content: space-around; margin-top: 10px;"> <div style="border: 1px solid black; padding: 5px; width: 45%;">                     A yellow precipitate appears; presence of lead cation  <math display="block">\text{Pb Cl}_2 + \text{K}_2\text{CrO}_4 \rightarrow 2\text{K Cl} + \text{Pb CrO}_4(\text{Yellow})</math> </div> <div style="border: 1px solid black; padding: 5px; width: 45%;">                     No yellow precipitate appears                      Absence of lead cation                 </div> </div> 2. To the second portion of the centrifugate, potassium iodide solution is added $\text{PbCl}_2 + 2 \text{KI} \rightarrow 2\text{KCl} + \text{PbI}_2$ (golden spangles) A yellow precipitate appears                      It is boiled with water to dissolve and cooled under tap water <div style="display: flex; justify-content: space-around; margin-top: 10px;"> <div style="border: 1px solid black; padding: 5px; width: 45%;">                     Golden spangles appear                      Presence of lead cation <math>\text{Pb Cl}_2 + 2 \text{K Cl} \rightarrow \text{Pb I}_2</math> </div> <div style="border: 1px solid black; padding: 5px; width: 45%;">                     No golden spangles appear                      Absence of lead cation                 </div> </div>

### ANALYSIS OF GROUP II CATIONS

The residue for group II is washed with water to remove alkali and the washings are discarded. Then a few drops of dilute $\text{HNO}_3$ are added and boiled. Then a few drops of dilute $\text{H}_2\text{SO}_4$ are added and centrifuged. $\text{CuS}$ or $\text{Cd S}$ or $\text{Bi}_2\text{S}_3 + \text{HNO}_3 \rightarrow$ Soluble nitrate of group II cations	
<b>NO RESIDUE</b> Absence of Mercury and Lead cations	<b>CENTRIFUGATE:</b> To the centrifugate, aqueous ammonia is added in drops to slight excess, heated and centrifuged $\text{Nitrates of group II cations} + \text{NH}_4 \text{OH} \rightarrow [\text{Cu}(\text{NH}_3)_4]^{2+} \text{ aqueous} + [\text{Cd}(\text{NH}_3)_4]^{2+} \text{ aqueous} + \text{Bi}(\text{OH})_3 \text{ residue.}$
<b>RESIDUE</b> It is dissolved by adding dilute $\text{HCl}$ and divided into two portions. $\text{Bi}(\text{OH})_3 + 3 \text{HCl} \rightarrow 3 \text{H}_2\text{O} + \text{BiCl}_3$	<b>CENTRIFUGATE</b> To the centrifugate, dilute $\text{H}_2\text{SO}_4$ is added followed by hydrogen sulphide gas. The precipitate is heated with a few drops of conc. $\text{HCl}$ and again $\text{H}_2\text{S}$ gas is passed and centrifuged $\text{Cu SO}_4 + \text{H}_2\text{S} \rightarrow \text{H}_2\text{SO}_4 + \text{Cu S (black precipitate)}$

	<p>To the first portion, a few drops of sodium hydroxide and sodium stannite solution are added.</p> $2 \text{Bi}^{3+} + 6 \text{OH}^- + 3 \text{Na}_2 \text{SnO}_2 \rightarrow 3 \text{Na}_2 \text{SnO}_3 + 3 \text{H}_2\text{O} + 2 \text{Bi (black ppt.)}$	<p><b>RESIDUE</b></p> <p>The black residue is dissolved in drops of dilute <math>\text{HNO}_3</math> and 3 props of dilute acetic acid are added. Then potassium ferrocyanide solution is added.</p> $\text{Cu S} + 2 \text{HNO}_3 \rightarrow 2 \text{H}^+ + \text{Cu(NO}_3)_2$ $\text{Cu}^{2+} + 2 \text{OH}^- \rightarrow \text{Cu(OH)}_2 \text{ (blue)}$ $\text{Cu}^{2+} + \text{K}_4 [\text{Fe(CN)}_6] \rightarrow \text{Cu}^{2+}[\text{Fe(Cn)}_6] \text{ (red)}$	<p><b>CENTRIFUGATE</b></p> <p>The centrifugate is diluted with water and hydrogen sulphide is passed</p> $\text{Cd}^{2+} + \text{H}_2\text{S} \rightarrow 2 \text{H}^+ + \text{Cd S (yellow precipitate)}$
	<p>Black and white precipitate appears</p> <p>Presence of Bismuth cation.</p>	<p>Reddish-brown Precipitate appears</p> <p>Presence of copper cation</p>	<p>Yellow precipitate appears</p> <p>Presence of cadmium cation</p>
	<p>To the second portion, a slight excess of water is added</p> $\text{Bi}^{3+} + 3 \text{H}_2\text{O} \rightarrow 6\text{H} + \text{Bi(OH)}_3 \text{ precipitate}$	<p>No precipitate appears</p> <p>Absence of copper cation</p>	<p>No precipitate appears</p> <p>Absence of cadmium cation</p>
	<p>White turbidity appears.</p> <p>Presence of Bismuth cation</p>		

### ANALYSIS OF GROUP III CATIONS

The residue for group III cations is washed with water and boiled with a little of sodium peroxide solid. It is stirred and boiled for some more time and centrifuged.  $Mn^{2+}$  or  $Fe^{2+}$  or  $Al^{3+} + Na_2O_2 \rightarrow Na_3 [Al(OH)_6]$  (in centrifugate)

<p><b>RESIDUE:</b></p> <p>The residue is dissolved dilute HCl and divided into two parts          To the first portion ammonium thiocyanate solution is added  <math>Fe(OH)_3 + 3 HCl \rightarrow Fe Cl_3 + 3 H_2O</math>  <math>Fe^{3+} + 3 NH_4 CNS \rightarrow 3NH_4Cl + Fe(CNS)</math> (bloodred colour)</p> <div style="display: flex; justify-content: space-between; margin-top: 10px;"> <div style="width: 45%;"> <p>Blood-red colour appears Presence of iron cation</p> </div> <div style="width: 45%; border: 1px solid black; padding: 2px; text-align: center;"> <p>No colour appears Absence of iron cation</p> </div> </div> <p style="margin-top: 10px;">To the second portion, potassium ferrocyanide solution is added  <math>FeCl_3 + K_4 [Fe(CN)_2] \rightarrow 3Cl^- + K Fe[Fe(CN)_6]</math> (Prussian blue)</p> <div style="display: flex; justify-content: space-between; margin-top: 10px;"> <div style="width: 45%; border: 1px solid black; padding: 2px;"> <p>Blue colour appears Presence of iron cation</p> </div> <div style="width: 45%; border: 1px solid black; padding: 2px; text-align: center;"> <p>No blue colour appears Absence of iron cation</p> </div> </div>	<p><b>CENTRIFUGATE:</b></p> <p>It is divided into two portions          The first portion is boiled with solid ammonium chloride  <math>Al^{3+} + 3 NH_4OH \rightarrow 3 NH_4^+ + Al(OH)_3</math> (gelatinous precipitate)</p> <p style="margin-top: 10px;">Gelatinous white precipitate appears Presence of aluminium cation</p> <p style="margin-top: 10px;">2. To the second portion, dilute HCl is added followed by aluminon reagent and aqueous ammonia  <math>Na_2 [Al(OH)_6] + 6HCl \rightarrow 3 Na Cl + 6H_2O + Al Cl_3</math>  <math>Al Cl_3 + 3 NH_4 OH \rightarrow 3 NH_4Cl + Al(OH)_3</math>          This <math>Al(OH)_3</math> absorbs the coloring matter and forms red-lake</p> <div style="display: flex; justify-content: space-between; margin-top: 10px;"> <div style="width: 45%; border: 1px solid black; padding: 2px;"> <p>Red-lake appears Presence of aluminum cation</p> </div> <div style="width: 45%; border: 1px solid black; padding: 2px; text-align: center;"> <p>No characteristic reaction Absence of aluminum cation</p> </div> </div>
---	--

### ANALYSIS OF GROUP IV CATIONS

The residue for group IV cations, is washed with water and digested with a little of dilute HCl acid and centrifuged Zn S and Mn S are soluble in dilute HCl and they are present in the centrifugate, leaving CoS and NiS residue		
<b>RESIDUE</b> The residue containing CoS and NiS is dissolved in least quantity of conc. HNO <sub>3</sub> . It is evaporated just to dryness. The final residue is dissolved in 1 ml of water and divided into 3 portions.	<b>CENTRIFUGATE</b> It is boiled to expel hydrogen sulphide gas and sodium hydroxide is added in excess. Then it is centrifuged. $\text{Zn Cl}_2 + 4 \text{ NaOH} \rightarrow \text{Na}_3 [\text{Zn}(\text{OH})_4] + 2 \text{ Na Cl}$ $\text{Mn Cl}_2 + 2 \text{ NaOH} \rightarrow \text{Mn} (\text{OH})_2 + 2\text{Na Cl}$	
To the first portion, a few drops of ammonium thiocyanate and amyl alcohol are added and shaken well. $\text{CoCl}_2 + 4 \text{ NH}_4 \text{ CNS} \rightarrow 2 \text{ NH}_4 \text{ Cl} + (\text{NH}_4)_2 [\text{Co}(\text{CNS})_4] \text{ [deep blue complex]}$ <div style="border: 1px solid black; padding: 2px; margin: 5px 0;">                     Blue layer appears                      Presence of cobalt cation                 </div> To the second portion, ammonium chloride, ammonium hydroxide and potassium ferricyanide solutions are added one by one. $\text{CoCl}_2 + \text{K}_3 [\text{Fe}(\text{CN})_6] \rightarrow \text{Co}(\text{CN})_2 \text{ [reddish-brwn precipitate]}$ <div style="border: 1px solid black; padding: 2px; margin: 5px 0;">                     Reddish brown precipitate appears                      Presence of cobalt cation                 </div> To the third portion, ethanol, DMG and aqueous ammonium hydroxide solution are added Rosy Read precipitate or colour is formed. Presence of nickel cation	<b>RESIDUE</b> It is dissolved in 1 ml of dilute HNO <sub>3</sub> and divided into two portions To the first portion sodium bismuthate is added stirred and centrifuged. <div style="border: 1px solid black; padding: 2px; margin: 5px 0;">                     Pink colour appears                      Presence of Mn cation                 </div> 2. To the second portion a pinch of Lead dioxide is added, boiled and allowed to stand <div style="border: 1px solid black; padding: 2px; margin: 5px 0;">                     Pink colour appears                      Presence of Mn cation                 </div>	<b>CENTRIFUGATE</b> It is divided into two portions To the first portion, hydrogen sulphide gas is passed $\text{Na}_3 [\text{Zn}(\text{OH})_4] + 2\text{H}_2\text{S} \rightarrow \text{Na}_2\text{S} + 4 \text{ H}_2\text{O} + \text{Zn S}$ Dirty white precipitate appears It is dissolved in conc. HNO <sub>3</sub> after removing H <sub>2</sub> S and cobalt nitrate solution is added. A strip of filter paper is dipped in the solution and burnt to an ash. <div style="border: 1px solid black; padding: 2px; margin: 5px 0;">                     Green ash appears .                      Presence of Zn cation                 </div> 2. To the second portion, a few drops of dilute HOAc are added followed by potassium ferrocyanide. $\text{Na}_2[\text{Zn}(\text{OH})_4] + 4 \text{ HOAc} \rightarrow \text{Zn} (\text{OAc})_2 + 2\text{Na OAc} + 4 \text{ H}_2\text{O}$ $2 \text{ Zn}(\text{OAc})_2 + \text{K}_4[\text{Fe}(\text{CN})_6] \rightarrow 4 \text{ K OAc} + \text{Zn}_2\text{Fe}(\text{CN})_6$ [dirty white precipitate ] Dirty white precipitate appears Presence of Zinc cation.

## ANALYSIS OF GROUP V CATIONS

The residue for group V cations, is washed with water and dissolved in minimum quantity of dilute acetic acid and 1 ml of potassium chromate is added for complete precipitation and centrifuged. [carbonates of barium, strontium and calcium] $MCO_3 + 2 HOAc \rightarrow M(OAc)_2 + CO_2 + H_2O$ $Ba(OAc)_2 + K_2 CrO_4 \rightarrow 2 K OAc + Ba CrO_4 \text{ (precipitated leaving the carbonates of Sr and Ca in centrifugate)}$			
<b>RESIDUE</b> Yellow residue is mixed with a drop of conc. HCl and made as a paste. The paste is introduced at the base of the outer oxidizing region of the non-luminous flame of burner.	<b>CENTRIFUGATE:</b> The centrifugate is neutralized with aqueous ammonia and ammonium carbonate is added and centrifuged. The residue is dissolved in dilute acetic acid and divided into two portions.		
<div style="border: 1px solid black; padding: 2px; margin-bottom: 10px;">                     Green Flame appears                      Presence of barium cation                 </div> <div style="border: 1px solid black; padding: 2px;">                     No Flame appears                      Absence of barium cation                 </div>	1. To the first portion, 1 ml of conc. $H_2SO_4$ are added. $Sr CO_3 + H_2SO_4 \rightarrow Sr SO_4$  <div style="border: 1px solid black; padding: 2px; margin: 10px auto; width: 80%;">                     White precipitate appears                      Presence of strontium cation                 </div>	2. To the second portion, saturated solution of ammonium sulphate is added and centrifuged. To the centrifugate, 1 ml of HOAc and ammonium oxalate solutions are added followed by aqueous ammonia $CaCO_3 + (NH_4)_2 SO_4 \rightarrow (NH_4) CO_3 + Ca SO_4$ $Ca SO_4 + 2 HOAc \rightarrow H_2SO_4 + Ca (OAc)_2$ $Ca(OAc)_2 + (NH_4)_2 C_2 O_4 \rightarrow 2NH_4 OAc + CaC_2 O_4$	
		White precipitate appears Presence of calcium cation	No white precipitate appears Absence of calcium cation
<b>ANALYSIS OF GROUP VI CATIONS</b>			
The centrifugate from group V is divided into two portions			
1. To the first portion ammonium chloride solution and aqueous ammonia are added followed by the addition of sodium hydrogen phosphate reagent. The sides of the test tube is scratched with a glass rod. $Mg^{2+} + Na_2 H PO_4 + NH_3 \rightarrow 2 Na^+ + Mg NH_4 PO_4 \text{ (white precipitate)}$			
<div style="border: 1px solid black; padding: 2px; margin-bottom: 10px;">                     White precipitate appears                      Presence of magnesium cation                 </div>		<div style="border: 1px solid black; padding: 2px;">                     No white precipitate appears                      Absence of magnesium cation                 </div>	
2. To the second portion a few drops of dilute HCl are added, stirred and a drop of Magneson reagent is added. Then a few drops of sodium hydroxide are added.			
<div style="border: 1px solid black; padding: 2px; margin-bottom: 10px;">                     White precipitate appears                      Presence of magnesium cation                 </div>		<div style="border: 1px solid black; padding: 2px;">                     No blue precipitate appears                      Absence of Magnesium cation                 </div>	



### TEST FOR AMMONIUM CATION

Since many ammonium salts are added during the course of analysis ammonium cation has to be tested in the unknown before proceeding to the systematic analysis for cations.

The given unknown mixture is washed with water and centrifuged.

To the centrifugate, a few drops of Nessler's reagent are added followed by the addition of sodium hydroxide solution

A brown precipitate appears  
Presence of ammonium cation

No brown precipitate appears  
Absence of ammonium cation

REPORT The given unknown mixture contains the following anions and cation

ANIONS: (I) Non interfering \_\_\_\_\_ CATIONS (1) \_\_\_\_\_ (Group \_\_\_\_\_)  
(II) Interfering \_\_\_\_\_ (2) \_\_\_\_\_ (Group \_\_\_\_\_)